











FUNDAMENTALS OF X-RAYS.....	2
X-RAY SETTINGS.....	3
RADIATION DOSE.....	5
DOSIMETRY	5
EFFECTS	6
DOSE REDUCTION.....	6
SOURCES OF RADIATION EXPOSURE	6
FILM VS DIGITAL IMAGING	7
DETECTOR CHARACTERISTICS.....	8
TROUBLESHOOTING.....	8
TYPES OF RADIOGRAPHS	8
SPECIALIZED VIEWS.....	9
TECHNIQUES.....	9
COMMON ERRORS.....	10
RADIOGRAPH INTERPRETATION	10
BENIGN VS MALIGNANT LESIONS.....	11
CARIES DETECTION.....	12
PERIODONTITIS DETECTION.....	12



Fundamentals of X-Rays

Position-Indicating Device (PID)	<p>Allow you to line up the tube head and distance from the patient</p> 										
Within the Tube Head	<p>A circuit heats up a filament -> Emits electrons from the Cathode to the Anode Electrons fire towards a Tungsten target -> X-Ray is produced when e^- hits the atoms</p> <p>*Copper rod is present to dissipate the heat</p>   <table border="1" data-bbox="391 905 1511 1108"> <tr> <td>Cathode (-)</td><td>Tungsten filament that produces electrons <ul style="list-style-type: none"> - Molybdenum focusing cup focuses the beam onto a small focal spot </td></tr> <tr> <td>Anode (+)</td><td>Tungsten target converts electrons to X-ray photons <ul style="list-style-type: none"> - Copper stem dissipates heat </td></tr> <tr> <td>Glass</td><td>Insulation</td></tr> <tr> <td>Aluminum</td><td>Filtration</td></tr> <tr> <td>Lead</td><td>Collimation</td></tr> </table>	Cathode (-)	Tungsten filament that produces electrons <ul style="list-style-type: none"> - Molybdenum focusing cup focuses the beam onto a small focal spot 	Anode (+)	Tungsten target converts electrons to X-ray photons <ul style="list-style-type: none"> - Copper stem dissipates heat 	Glass	Insulation	Aluminum	Filtration	Lead	Collimation
Cathode (-)	Tungsten filament that produces electrons <ul style="list-style-type: none"> - Molybdenum focusing cup focuses the beam onto a small focal spot 										
Anode (+)	Tungsten target converts electrons to X-ray photons <ul style="list-style-type: none"> - Copper stem dissipates heat 										
Glass	Insulation										
Aluminum	Filtration										
Lead	Collimation										
What is an X-Ray?	<p>High Energy and High Frequency waves</p> <ul style="list-style-type: none"> - Between UV rays and Gamma Rays on the spectrum 										
Attenuation	<p>= How the X-Ray beam weakens as it passes through matter on its way to the receptor</p> <ul style="list-style-type: none"> - Density of the material it is passes through AND the energy of the photon (kVp) determines how many photons will actually make it to the receptor 										
Filtration	<p>= Involves Aluminum</p> <ul style="list-style-type: none"> - Removes the lower energy photos from the bean to \downarrow patient exposure (and \uparrow image quality) - Conceptually similar to beam hardening 										
Collimation	<p>= Involves Lead</p> <ul style="list-style-type: none"> - \downarrow Scatter radiation by \downarrow diameter of the beam - Rectangular is the best method to \downarrow patient exposure and also \uparrow image quality (by \downarrow scatter radiation)  										



2 Types of Ionizing Radiation

Electromagnetic	= Movement of energy as a combination of electrical and magnetic fields - Shorter Wavelength = ↑ Energy Gamma > X-Ray > UV > Violet-Red > Infrared > Microwave > Radio	
Particulate	= Atomic nuclei or subatomic particles moving at high velocity - Alpha (α) or Beta (β) particles from Radioactive decay - α is stopped by paper (pussy) - β is stopped by aluminum (this is why we use filtration)	


2 Types of X-Ray Production







Bremsstrahlung	= 1 ^o source of X-Ray photons - Electron strikes a target and brakes -> Gives off energy - Generates a continuous spectrum of energy	
Characteristic	= 2 ^o source of x-ray photos - Electron accidentally knocks into another electron and drops into a lower energy orbital - Emits a photon of specific energy when it moves orbitals	

X-Ray Settings

Definitions		
Intensity	= Quantity of electrons and photons Effects: Density - Darkness of the image	
Energy	= Quality of electrons and Energy of photos Effects: Contrast - Difference among gray values	

Variables




Exposure Time (s)	= Length of time that high voltage is applied to the X-Ray tube and the time that the current flows and produces X-Rays - Changes the NUMBER of photons - The most frequently altered setting If: Too long -> Too dark, overexposed Too short -> noisy and light, underexposed Affects: - Intensity (Density)	
--------------------------	--	---

Tube Current (mA)	<p>= Flow of e^- from filament to the anode and back to the filament</p> <ul style="list-style-type: none"> - Not really adjustable easily <p>If:</p> <ul style="list-style-type: none"> - Too much -> Too dark, underexposed - Too little -> Noisy, underexposed <p>Affects:</p> <ul style="list-style-type: none"> - Intensity (Density) 	
Tube Potential (kVp)	<p>= acceleration of e^- from cathode to anode</p> <ul style="list-style-type: none"> - Affects both the number and energy of the photons <p>If:</p> <ul style="list-style-type: none"> - Too high -> Too gray (not enough contrast), mostly Compton Scattering - Too Low -> Too light (very high contrast), Mostly photoelectric absorption <p>Affects:</p> <ul style="list-style-type: none"> - Density and Intensity 	
Distance	<p>*Operator should be at least 6ft away*</p>  <p>Ideal = ↑ Source Object distance and ↓ Object-Image Distance</p> <ul style="list-style-type: none"> - ↑ SOD = ↓ intensity (density), ↑ Quality - ↓ OID = ↑ Quality, ↓ Magnification <p>Geometry</p> <ul style="list-style-type: none"> - Umbra = Shadow behind an image - Penumbra = Side Shadow (blurriness) <ul style="list-style-type: none"> - Size depends on: Focal Spot Size, Source-to-Object distance, Object to image distance - Ideal: ↓ FSS, ↓ OID, and ↑ SOD = Chris P image 	  

Everything in one Table ->

	Intensity	Energy
↑ seconds	↑ density	
↑ mA	↑ density	
↑ kVp	↑ density	↓ contrast
↑ filter	↓ density	↑ contrast
↑ distance	↓ density	

Radiation Dose



Coherent Scattering	<p>= Incident photon contacts an outer electron</p> <ul style="list-style-type: none"> - No energy loss, just a change in direction - Only about 8% of interactions (pretty minor) <p>Effect:</p> <ul style="list-style-type: none"> - ↓ Contrast 	
Photoelectric Absorption	<p>= Incident photon contacts and ejects an inner electron and forms an ion pair</p> <ul style="list-style-type: none"> - If your kVp is too ↓ ↓ you will get ↑ Photoelectric Absorption - About 30% of interactions <p>Effect:</p> <ul style="list-style-type: none"> - ↑ Contrast 	
Compton Scattering	<p>= Incident photon contacts an outer electron and forms an ion pair</p> <ul style="list-style-type: none"> - If kVp is ↑ ↑ you get ↑ Compton Scattering - About 62% of interactions in a dental X-ray beam <p>Effects:</p> <ul style="list-style-type: none"> - ↓ Contrast 	

Dosimetry

Term	Unit of Measure	What It Is
<u>Exposure</u>	Roentgen 1 C/kg = 3.88 x 10³ R	Energy produced by x-ray tube
Absorbed Dose	Gray 1 Gray = 100 rads	Energy absorbed in tissue
Equivalent Dose	Sievert 1 Sievert = 100 rems	Energy absorbed in tissue multiplied by radiation weighting factor
Effective Dose	Sievert 1 Sievert = 100 rems	Energy absorbed in tissue multiplied by tissue weighting factor

→ Different tissues are affected differently (↑/↓)

Effects

Deterministic	<p>= Hair loss, Cataracts, Skin Damage, Oral Mucositis</p> <p>Threshold Dose needs to be reached -> ↓ than this dose = no effects</p> <ul style="list-style-type: none"> - 0.1Gy = In-utero birth defects - 0.5 Gy = Cataracts - 3.0 Gy = Radiation burns <p>**Dental X-Rays are Micro Grays, so they are FAR FAR away from causing Deterministic Effects**</p>	
Stochastic Effects	<p>= Cancer, Leukemia, Heritable Effects</p> <ul style="list-style-type: none"> - Accumulates over years of doses. This is why we are conservative <p>Linear, no-threshold model. Dose builds up over years</p> <ul style="list-style-type: none"> - 0.1Gy -> Anything below this, we can only assume the effects are dependent on dose (Extrapolating) - 0.0017% risk ↑ by 1Gy 	
Radiation Chemistry		
Direct	<p>= Direct alteration of biologic molecules (Proteins, DNA etc)</p> <ul style="list-style-type: none"> - Accounts for 1/3 of biologic effects - Effects caused by α and β particles 	
Indirect	<p>= Ionizing radiation converts H₂O into free radicals which then alter biologic molecules (proteins, DNA etc)</p> <ul style="list-style-type: none"> - Account for 2/3 of biologic effects <p>**Cells that are mitotically active are more radiosensitive**</p> <ul style="list-style-type: none"> - Most Sensitive: Hematopoietic Cells, Epithelial Cells, Sperm Cells - Least Sensitive: Nerve cells, Muscle cells 	

Dose Reduction

- Occupational Exposure limit = 50mSv per year
- ALARA = As Low As Reasonably Achievable
- ALADA = As Low As Diagnostically Acceptable
- Use E/F speed films, or digital imaging
- >18cm source to object distance to ↓ exposure and ↑ clarity
- Rectangular Collimation


Sources of Radiation Exposure


- Background radiation (just surviving on earth) = 3.1mSv per year -> mostly from Radon
- Man-made radiation = 3.1mSv per year -> Mostly medical imaging, building materials, etc
- Average Annual exposure = 6.2 mSv

****Dental X-Rays are MICRO Sv...SOOOOO much ↓ ↓ ↓ than what we are exposed to on a daily basis**

Technique	Effective dose in μSv	Dose as multiple of average † panoramic dose	Dose per capital background	Probability of x in a million fatal cancer†
FMX with PSP or F Speed film w Rectangular Collimation	35	2.2	4.3	2
FMX with PSP or F Speed film w Round Cone	171	10.7	21	9
FMX with D Speed film w round cone†	388	24.1	47	21
Panoramic – CCD†	16.1	1.0	2	0.9
Lateral Cephalometric - PSP	5.6	0.3	0.7	0.3
CB Mercuray- "Panoramic" FOV	560	35	68	30.8
Galileos – (default exposure)	70	4.3	9	3.9
Somaton 64 MDCT	860	53.4	105	47.3

Film vs Digital Imaging





Film	Digital
<ul style="list-style-type: none"> - Requires Chemicals to process - Requires time to develop - Superior image quality - More radiation dose to the patient <p>Composition:</p> <ul style="list-style-type: none"> - <u>Base</u> = Flexible Plastic - <u>Emulsion</u> = Silver halide crystals in gelatin material (sensitive to light) -> this is where the magic happens - <u>Intensifying Screens</u> = Coated w/ fluorescent phosphor to ↓ amount of exposure needed (converts photon energy to light), but also ↓ resolution  <p>Speed:</p> <ul style="list-style-type: none"> - $A < B < C < D < E < F$ - <u>Determined by several factors:</u> <ul style="list-style-type: none"> - Larger crystals = faster films - Double emulsion = Faster film - Radiosensitive dyes in emulsion = Faster film <p>Imaging:</p> <ul style="list-style-type: none"> - X-Ray photos chemically change the silver halide crystals into neutral silver atoms in the emulsion -> Creates latent image <p>Chemical Processing</p> <ol style="list-style-type: none"> 1. Developing 2. Fixing <ul style="list-style-type: none"> - Wash away unexposed silver halide crystals 3. Washing <ul style="list-style-type: none"> - Wash away residual chemicals 4. Drying 	<ul style="list-style-type: none"> - No chemicals needed - Instant viewing - Image enhancement w/ computer software - ↓ Radiation Dose to patient <p>Ex:</p> <ul style="list-style-type: none"> - PSP - CCD/CMOS <p>PSP:</p> <ul style="list-style-type: none"> - Barium Fluorohalide plates capture and store X-ray energy from dental exposure - Need scanner to scan the image into the computer <p>CCD/CMOS:</p> <ul style="list-style-type: none"> - Silicon sensor chip captures x-ray and rapidly displays image on the monitor

Chemical Processing	
Developing	<p>= Converts exposed crystals into metallic silver grains that appear dark on the radiograph</p> <ul style="list-style-type: none"> - Phenidone = 1st e⁻ donor that reduces silver ions to metallic silver at latent image site - Hydroquinone = Provides electron to reduce the oxidized phenidone back to its original active state (so it can continue to work) <p>If final image is light: Developer is probably old</p> 
Fixing	<p>= Wash away unexposed silver halide crystals</p> <ul style="list-style-type: none"> - Ammonium Thiosulphate = Cleaning agent, removes undeveloped silver halide crystals - Aluminum Salts = Tanning agent, hardens and preserves the emulsion - Acetic Acid = Acidifier, maintains the fixer's level of acidity and neutralizes the developer - Sodium Sulfite = Preservative, Prolongs shelf life of the solution - Water = Solvent, Dissolves the other ingredients
Washing	= Wash away residual chemicals
Drying	

Detector Characteristics







Contrast Resolution	= Ability to distinguish shades of Gray Film > Digital
Spatial Resolution	= Ability to distinguish 2 close points Film > CCD/CMOS > PSP
Detector Latitude	= Exposure range providing useful image PSP > CCD/CMOS > Film
Detector Sensitivity	= Dose required to achieve standard gray level CCD/CMOS is ½ the speed of n F Speed Film -> WAY less radiation

Troubleshooting

Underexposed Image	<ul style="list-style-type: none"> - Inadequate Exposure Time - Inadequate Development Time - Old Development Solution - Temp too low 
Overexposed Image	<ul style="list-style-type: none"> - Exposed for too long - Exposed to light 
Creases	
Static Electricity	
Film is placed backwards	= Herringbone/Tire Track  Effect




Types of Radiographs

Intraoral/Direct	Extraoral/Indirect
<ul style="list-style-type: none"> - Periapical (PA) - Bitewing (BW) - Occlusal 	<ul style="list-style-type: none"> - Panoramic (Pan) - Cephalometric (Ceph) - Cone Beam Computed Tomography (CBCT)



Intraoral	
Periapical	= Used to capture the roots and apices of teeth <ul style="list-style-type: none"> - Manifestations of pulpal or periodontal disease = Used also to capture anterior interproximal decay  
Bitewing	*Does not usually show root apices = Used for posterior interproximal caries detection and bone levels <ul style="list-style-type: none"> - Angulation is more parallel for BW's so bone is more accurate  
Occlusal	Receptor lies flat on the occlusal plane <ul style="list-style-type: none"> - Captures from Canine to Canine = Used for further assessments of: Trauma, Impacted Teeth, Supernumerary  

Extraoral	
Panoramic Cephalometric	<p>Useful for screening, pathology, and 3rd molar location</p> <p><u>Lateral</u></p> <ul style="list-style-type: none"> - Relationship of teeth and jaws in relation to the cranium - Treatment changes - Growth <p><u>Posterior-Anterior (PA)</u></p> <ul style="list-style-type: none"> - Transverse Asymmetry - Not really used anymore, mostly CBCT now
Cone Beam Computed Tomography (CBCT)	<p><u>3D radiographic image:</u></p> <ul style="list-style-type: none"> - Axial View (top bottom) - Sagittal (Side-side view) - Coronal (Front Back view) - 3D Volumetric render <p><u>Used for:</u></p> <ul style="list-style-type: none"> - Implant Planning - Endodontics (strange anatomy) - Orthodontics - TMJ (Hard tissue only) -> MRI is used for soft tissue - Pathology









Specialized Views

Waters View	<p>PA Ceph at an angle (45° to the orbito medial line)</p> <p><u>Best for:</u></p> <ul style="list-style-type: none"> - Paranasal Sinus - Midface - Orbits 	
Towne's View	<p>PA Ceph at an angle again (30° from the orbito-medial line)</p> <ul style="list-style-type: none"> - Angled directly at the condyles. Eliminates the superimposition of the Mastoid and Zygoma over the condyles <p><u>Best for:</u></p> <ul style="list-style-type: none"> - Condyle 	
Submentovertex View	<p>Base projection of the skull</p> <p><u>Best for:</u></p> <ul style="list-style-type: none"> - Basilar skull fractures - Zygomatic fractures 	

Techniques

Bisecting Angle Technique	<p>= Central ray of the X-Ray beam is aimed perpendicular to the imaginary bisecting line between the long axis of the tooth and the long axis of the receptor</p> <ul style="list-style-type: none"> - This creates 2 equal triangles...which means there is minimal magnification! 	
Paralleling Technique	<p>= Receptor is placed parallel to the long axis of the tooth</p> <ul style="list-style-type: none"> - Central ray of the x-ray beam is aimed perpendicular to the long axis of both the tooth and receptor <p>*Less distortion and ↑ quality image*</p> <ul style="list-style-type: none"> - Problem is that the Object film distance is ↑ (less comfy and more magnification) 	






Common Errors















Elongation	= Most common error - Could be angulation issue or bending of the film	
Cone Cut	x-Ray beam and receptor were not lined up properly	
Underexposed	Grainy or too light image - Exposure time is too short - ↓ mA - ↓ kVp	
Overexposed	Image too dark - ↑ mA - ↑ kVp	
Double Exposure	2 Images exposed on the same PSP plate	
Pans		
Motion	Wavy and Irregular borders - If patient swallows or moves	
Chin Down	= Big Smile	
Chin Up	= Frown	
Ghosting	RO Artifact is projected to the other side of the image	

Radiograph Interpretation

- Know what normal is
- Compare R and L sides (especially with Pan)
- Categorize the disease or abnormality
- Start globally -> then go locally


Descriptions

Radiolucent Descriptors					
Corticated Unilocular	One Compartment Radiopaque Border		Non-Corticated Unilocular	One Compartment No clear border Diffuse edges	
Multilocular	Multiple compartments separated by thin septae		Multifocal Confluent	Multiple points of origin converging onto each other	
Moth-Eaten	Irregular, ragged edges Can be generalized or localized				

Radiopaque Descriptors			
Focal Opacity 	Single Site Homogenous = Same density Heterogenous = Varying density 	Target Lesion 	Radiopaque Center RL Band around the RO center Corticated border around the RL halo 
Multifocal Confluent 	Multiple Points of origin converging on each other 	Irregular 	Ill defined Edges 
Ground Glass 	Fine granular/orange peel appearance 	Mixed Density 	Both RO and RL May or may not be corticated 
Soft Tissue Opacity 	Calcification embedded within the soft tissue 		

Soft Tissue RO in a Pan






Benign vs Malignant Lesions



	Benign	Malignant
Margins	Well-defined (narrow zone of transition), smooth, regular, corticated	Ill-defined (wide zone of transition), ragged, moth-eaten
Shape	Round, oval	Irregular
Internal Architecture	Tumors Multilocular, more likely to resorb roots Cysts Unilocular, hydrostatic, tend to be more corticated	Usually radiolucent
Location	Coronal or above mandibular canal	Ramus or posterior body of mandible
Effect on Cortical Bone	Expansion, thinning, displacement, may erode if aggressive	Erosion, destruction
Effect on Maxillary Sinus	Displacement	Erosion, destruction
Effect on Mandibular Canal	Displacement, no neurosensory defects	Invasion and destruction, anesthesia or paresthesia
Effect on Tooth Position	Displaced teeth, may prevent eruption	Floating teeth
Effect on Tooth Roots	Horizontal resorption	Vertical resorption resulting in spiked roots, or no resorption at all (floating teeth)
Effect on PDL Space and Lamina Dura	N/A	Asymmetric widening of PDL space and loss of lamina dura



Caries Detection

Tooth needs 30-40% demineralization before it will appear on a radiograph

Interproximal	Small triangle at or below the contact point		
Occlusal	Subtle RL beneath the fissure - Hard to trace		
Buccal	Small circle in the buccal or lingual pit		
Recurrent	Gingival to the restoration		
Root	Hemispherical at or below the CEJ - Important to distinguish from cervical burnout		

Periodontitis Detection

- Use BW Radiographs, the Bone levels are more accurate -> Except for Anteriors where we can't really use BW's
- Alveolar crest is fuzzy =
- Normal bone levels: 1-2mm apical to the CEJ

Mild Bone Loss	Upper 3 rd of root		
Moderate BL	Middle 3 rd of root		
Severe BL	Lower 3 rd of root	